3. EXPERIMENTAL RESULTS

3.1. Specific heat of solid helium

The specific heat of solid 4 He was measured at four different molar volumes, and of solid 3 He at ten different molar volumes. The results are shown in figures 5 and 6. Debye temperatures θ_{D} were calculated from the experimental points and are shown in figures 7 and 8. Included in figure 7 are earlier measurements by Keesom & Keesom (1936).




FIGURE 7. The Debye temperature of solid ⁴He as a function of temperature at different molar volumes.

The heat capacity of the helium samples contributed from 30 to 40 % of the total measured heat capacity near the melting point. This ratio became progressively worse with falling temperature and varied from 5 to 25 % at 3 °K. The accuracy is therefore greatly reduced below about 5 °K, especially for samples of low molar volume. Because of this we feel that